Calvian® / Calvian endo-pen®

Bipolare Koagulationszangen

Bipolar Coagulation Forceps

Pincas de Coagulation Bipolaires

Calvian® / Calvian endo-pen®

Pinzas de Coagulación Bipolares

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Bipolare Coagulationtangen

DE EN FR

ES IT NL

Deutsch

DE

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STERRAD® is a trademark of Johnson & Johnson, Inc.

English

EN

Product / User/ Disposal:

Electrosurgery accessories may only be used and disposed of by qualified medical personnel! Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.

Not sterile. Must be cleaned and sterilized prior to initial and each subsequent use.

Intended use:

Bipolar coagulation of soft tissue. Depending on the instrument's scope of functionality, also for suctioning fluids during surgical procedures.

Service life:

A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:

- Inspect the product for cleanliness, mechanical function, and intact insulation before each use. We recommend checking the insulation with a suitable test device. Only use sterilized products that are in flawless condition! A certain discoloration of the non-stick instrument tip is normal and harmless. To connect the instrument and cables, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

During use:

- Always work with the lowest output setting for the desired surgical effect. Maximum allowable voltage 500 Vp. Regularly wipe blood and tissue residue from the tips. Instrument tips may cause injuries! After application, instrument tips may be so hot they cause burns! Never lay down instruments on or in the immediate vicinity of the patient! Not for use in the presence of combustible or explosive substances!

Reconditioning:

General information:

Observe national guidelines and regulations! Disconnect the instrument from the cable! Do not allow blood and tissue residues dry up! Remove blood and tissue residues with a soft cloth or brush! Do not use sharp tools / scouring agents! Do not immerse in hydrogen peroxide (H2O2)!

Manual cleaning and disinfection:

The instrument must always be reconditioned by machine – no manual cleaning! According to recommendations of the DGSV ("Deutsche Gesellschaft für Sterilgutversorgung", German Association for Sterile Goods Supply), the Calvian® and Calvian endo-pen® are assigned to risk group B*. Machine cleaning is required for such products.

This classification was assigned according to the DGSV 2013 flow chart for the classification of medical devices based on the recommendations of the KRINKO ("Krankenhaushygiene und Infektionsprävention Kommission", Hospital Hygiene and Infection Prevention Commission)/BfArM ("Bundesinstitut für Arzneimittel und Medizinprodukte", Federal Institute for Drugs and Medical Devices), Bundesgesundheitsblatt (Federal Health Bulletin) 2012; 55:1244-1310

Manual pre-cleaning:

Manual pre-cleaning is an important prerequisite for successful machine cleaning and therefore part of the overall reconditioning process!

Table with 2 columns: Calvian® (1) and Calvian endo-pen® (2). Rows describe cleaning steps like 'Promptly after use', 'Clean the instrument (especially the instrument tip)', 'Thoroughly flush both lumina (channels)', and 'Rinse the lumen (channel)'.

Table with 2 columns: Calvian® (1) and Calvian endo-pen® (2). Row: Clean the instrument in an ultrasound bath: 40 °C, 15 minutes, mildly alkaline cleaner with a concentration of 0.5 % deconex® 28 ALKA ONE-x (Borer Chemie).

Proof of fundamental suitability of the instruments for manual pre-cleaning was provided by an independent, accredited test laboratory, test report no. 2105011203 dated 23.03.2005.

Proof of fundamental suitability of the instruments for manual pre-cleaning was provided by an independent, accredited test laboratory, test report no. 30514 dated 02.02.2015 / test report no. 03719 dated 26.03.2019.

Machine cleaning and disinfection:

In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883). The CDD must have connections for flushing the instruments.

Process:

Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported. The Calvian endo-pen® should be placed into the storage trays included or available as accessories (REF 701778-01, 701778-02, 701778-05 and 701778-10). Calvian endo-pen®, bayonet (REF 701778-01, 701778-02). The lumina of the instruments must be connected to the flushing connection of the CDD using the Luer-Lock connections provided.

Table with 3 columns: Program steps, Calvian® (1), and Calvian endo-pen® (2). Rows include 'Pre-rinse with cold tap water', 'Clean with 0.5 % deconex® 28 ALKA ONE-x at 70 °C', 'Neutralize with warm tap water', etc.

Proof of fundamental suitability of the instruments for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory using the Vario TD / Miele G7735 CD disinfectant (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent deconex® 28 ALKA ONE-x with a concentration of 0.5 % (Borer Chemie AG, Zuchwil, Switzerland), test report no. 2105011203 dated 23.03.2005.

Proof of fundamental suitability of the instruments for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory using the Miele G7836 CD disinfectant (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent deconex® 28 ALKA ONE-x with a concentration of 0.5 % (Borer Chemie AG, Zuchwil, Switzerland), test report no. 30514 dated 02.02.2015.

Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instrument and are permitted according to the A0 concept, for instance thermal disinfection at 90 °C, 5 minutes, comparable A0-value->3000. If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Please contact your supplier or hygiene officer in case of questions.

Inspection:

Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:

None

Packaging:

Package cleaned and disinfected storage trays in disposable sterilization packaging (single or double packaging), or place the instrument or tray with the cleaned and disinfected instruments into a suitable sterilization container meeting the following requirements:

- EN ISO/ANSI AAMI ISO 11607 • Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability) • Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:

Products must be cleaned and disinfected prior to sterilization. Use only the sterilization process described in the following:

- Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Table with 3 columns: Program parameters, Calvian® (1), and Calvian endo-pen® (2). Rows include 'Process', 'Sterilization temperature', 'Sterilization time', 'Max. sterilization temperature', and 'Drying time'.

Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an independent, accredited test laboratory using the steam sterilizer Autoclave 6-6-6 Selectomat HP, test report no. 2105021003 dated 23.03.2005.

Proof of fundamental suitability of the instruments for effective steam sterilization was provided by an independent, accredited test laboratory using the steam sterilizer Autoclave 6-6-6 Selectomat 1HR, drying time min. 1 minute in a half cycle – here specification of a commonly used minimum drying time (e.g. for Melag Euroklav 23V-S 13-15 minutes "quick program not packaged"), test report no. 13514 dated 05.08.2014.

Do not sterilize in hot air!

Do not sterilize in STERRAD®!

In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Storage / Transportation:

Store in a cool, dry place. Protect from sunlight. Store and transport in secure containers / packaging. For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:

Serious incidents that occur in relation to the product must be reported to the manufacturer and the responsible public authority in the member state in which the user and/or patient resides.

The product may only be repaired by the manufacturer or an agent expressly authorized by the manufacturer. Otherwise the warranty and any possible additional liability claims against the manufacturer are voided.

Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.

Changes reserved. Current version available at www.sutter-med.de.

STERRAD® is a trademark of Johnson & Johnson, Inc.

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