

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Description:	DEXON™ II coated, braided synthetic absorbable sutures are composed of the homopolymer of glycolic acid and coated with polycaprolate, a copolymer of caprolactone and glycolide.
Characteristics:	The advanced manufacturing process of the molecule of polyglycolic acid and the coating with polycaprolate give the suture: - Excellent strength over the critical wound healing period - Excellent flow through tissue - Uniform diameter - Predictable absorption profile
Indications:	DEXON™ II sutures are indicated for use in soft tissue approximations or ligation including use in ophthalmic procedures, but not for use in cardiovascular tissue and/or neural tissue.
Actions:	DEXON™ II sutures elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of DEXON™ II sutures occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and metabolized by the body.
Tensile Strength:	Subcutaneous tissue implantation studies in rats indicate DEXON™ II sutures retain approximately 65% of the original tensile strength at two weeks and approximately 35% at three weeks post-implant.
Absorption:	Absorption begins as a loss of tensile strength without appreciable loss of mass. Absorption of DEXON™ II sutures is essentially complete between 60 and 90 days.
Sterilization:	Ethylene oxide.
Packaging:	DEXON™ II sutures are available sterile and in various lengths, diameters (USP/EP) and quantities with or without surgical needles. DEXON™ II sutures are undyed or dyed with D&C Green No. 6, D&C Violet No. 2, or bicolored.

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Non Needed								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
Pre-Cut Suture	4-0	1.5	12x45	Bi-Color	36			9805-31
	4-0	1.5	12x45	Bi-Color	24	UNI-SPENSE™		9809-32
	4-0	1.5	3x45	Bi-Color	36			9803-31
	4-0	1.5	3x45	Bi-Color	24	UNI-SPENSE™		9017-32
	4-0	1.5	2x75	Bi-Color	36		•	9802-31
	3-0	2	12x45	Bi-Color	36			9805-41
	3-0	2	12x45	Bi-Color	24	UNI-SPENSE™		9809-42
	3-0	2	3x45	Bi-Color	36			9803-41
	3-0	2	3x45	Bi-Color	24	UNI-SPENSE™		9017-42
	3-0	2	6x45	Bi-Color	24	UNI-SPENSE™		9004-42
	3-0	2	2x75	Bi-Color	36			9802-41
	3-0	2	7x75	Bi-Color	36			9806-41
	3-0	2	7x75	Bi-Color	24	UNI-SPENSE™		9877-42
	2-0	3	12x45	Bi-Color	36			9805-51
	2-0	3	12x45	Bi-Color	24	UNI-SPENSE™		9809-52
	2-0	3	3x45	Bi-Color	36			9803-51
	2-0	3	3x45	Bi-Color	24	UNI-SPENSE™		9017-52
	2-0	3	6x45	Bi-Color	24	UNI-SPENSE™		9004-52
	2-0	3	2x75	Bi-Color	36			9802-51
	2-0	3	7x75	Bi-Color	36			9806-51
	2-0	3	7x75	Bi-Color	24	UNI-SPENSE™		9877-52
	0	3.5	12x45	Bi-Color	36			9805-61
	0	3.5	12x45	Bi-Color	24	UNI-SPENSE™		9809-62
	0	3.5	3x45	Bi-Color	36			9803-61
	0	3.5	3x45	Bi-Color	24	UNI-SPENSE™		9017-62
	0	3.5	6x45	Bi-Color	24	UNI-SPENSE™		9004-62
	0	3.5	2x75	Bi-Color	36			9802-61
	0	3.5	7x75	Bi-Color	24	UNI-SPENSE™		9877-62
	1	4	12x45	Bi-Color	24	UNI-SPENSE™		9809-72
	1	4	12x45	Bi-Color	12			9805-73
	1	4	3x45	Bi-Color	36			9803-71
	1	4	3x45	Bi-Color	24	UNI-SPENSE™		9017-72
	1	4	6x45	Bi-Color	24	UNI-SPENSE™		9004-72
	1	4	2x75	Bi-Color	36			9802-71
	1	4	7x75	Bi-Color	24	UNI-SPENSE™		9877-72
	2	5	12x45	Bi-Color	24	UNI-SPENSE™		9809-82
	2	5	12x45	Bi-Color	12			9805-83
	2	5	3x45	Bi-Color	36			9803-81
	2	5	3x45	Bi-Color	24	UNI-SPENSE™		9017-82


Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Non Needled							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
Pre-Cut Suture	2	5	7x75	Bi-Color	24	UNI-SPENSE™	9877-82
Standard Length	4-0	1.5	150	Bi-Color	36		9801-31
	3-0	2	150	Bi-Color	36		9801-41
	2-0	3	150	Bi-Color	36		9801-51
	0	3.5	150	Bi-Color	36		9801-61
	1	4	150	Bi-Color	36		9801-71
	2	5	150	Bi-Color	36		9801-81
Ligating Reel	4-0	1.5	150	Bi-Color	12		9804-33
	4-0	1.5	250	Bi-Color	72		9808-37
	4-0	1.5	250	Bi-Color	12		9808-33
	3-0	2	150	Beige	24		9008-42
	3-0	2	150	Bi-Color	12		9804-43
	3-0	2	250	Bi-Color	72		9808-47
	3-0	2	250	Bi-Color	12		9808-43
	3-0	2	360	Bi-Color	12		9807-43
	2-0	3	150	Beige	24		9008-52
	2-0	3	150	Bi-Color	12		9804-53
	2-0	3	250	Bi-Color	72		9808-57
	2-0	3	250	Bi-Color	12		9808-53
	2-0	3	360	Bi-Color	12		9807-53
	0	3.5	150	Beige	24		9008-62
	0	3.5	150	Bi-Color	12		9804-63
	0	3.5	250	Bi-Color	72		9808-67
	0	3.5	250	Bi-Color	12		9808-63
	0	3.5	360	Bi-Color	12		9807-63
	1	4	150	Bi-Color	12		9804-73
	1	4	250	Bi-Color	72		9808-77
	1	4	250	Bi-Color	12		9808-73
	2	5	250	Bi-Color	12		9808-83




Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
3/8 CIRCLE 12 mm 	6-0	0.7	45	Beige	36		9538-11	
	5-0	1	45	Beige	36		9538-21	
	4-0	1.5	45	Beige	36		9538-31	
3/8 CIRCLE 16 mm 	5-0	1	75	Beige	36		9602-21	
	4-0	1.5	75	Beige	36		9602-31	
	3-0	2	75	Beige	36		9602-41	
3/8 CIRCLE 19 mm 	5-0	1	45	Beige	36		9608-21	
	5-0	1	75	Beige	36		9612-21	
	4-0	1.5	45	Beige	36		9608-31	
	4-0	1.5	75	Beige	36		9612-31	
	3-0	2	45	Beige	36		9608-41	
	3-0	2	75	Beige	36		9612-41	
3/8 CIRCLE 24 mm 	4-0	1.5	75	Beige	36		9618-31	
	4-0	1.5	75	Bi-Color	36		9619-31	
	3-0	2	75	Beige	36		9618-41	
	3-0	2	75	Bi-Color	36		9619-41	
	2-0	3	75	Beige	36		9618-51	
	2-0	3	75	Bi-Color	36		9619-51	
	2-0	3	75	Violet	36		9641-57	
3/8 CIRCLE 24 mm 	2-0	3	150	Bi-Color	36		9839-51	
1/2 CIRCLE 37 mm 								


Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
3/8 CIRCLE 30 mm 	3-0	2	75	Beige	36			9626-41
	2-0	3	75	Beige	36			9626-51
	2-0	3	75	Bi-Color	36			9627-51
	0	3.5	75	Beige	36			9626-61
	0	3.5	75	Bi-Color	36			9627-61
	1	4	75	Beige	36			9626-71
3/8 CIRCLE 39 mm 	2-0	3	75	Beige	36			9628-51
1/2 CIRCLE 20 mm 	4-0	1.5	75	Bi-Color	36			9860-31
1/2 CIRCLE 24 mm 	3-0	2	45	Beige	36			9735-41
	3-0	2	75	Beige	36			9740-41
	3-0	2	75	Bi-Color	36			9875-41
	2-0	3	75	Beige	36			9740-51
	2-0	3	75	Bi-Color	36			9875-51
1/2 CIRCLE 26 mm 	3-0	2	75	Bi-Color	36			9861-41
	2-0	3	5x45	Beige	36	D-TACH™		9746-51
	2-0	3	75	Bi-Color	36			9861-51
	0	3.5	75	Bi-Color	36			9861-61



Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 26 mm H 	2-0	3	5x45	Beige	36	D-TACH™		9520-51
	0	3.5	75	Bi-Color	36			9862-61
	1	4	75	Bi-Color	36			9862-71
1/2 CIRCLE 30 mm 	3-0	2	75	Bi-Color	36			9584-41
	2-0	3	75	Bi-Color	36			9584-51
	0	3.5	75	Bi-Color	36			9584-61
	1	4	75	Bi-Color	36			9584-71
1/2 CIRCLE 37 mm 	3-0	2	75	Bi-Color	36			9863-41
	2-0	3	75	Bi-Color	36			9863-51
	2-0	3	90	Bi-Color	36			9864-51
	0	3.5	75	Bi-Color	36			9863-61
	0	3.5	90	Bi-Color	36			9864-61
	1	4	75	Bi-Color	36			9863-71
	1	4	90	Bi-Color	36			9864-71
	2	5	75	Bi-Color	36			9863-81
1/2 CIRCLE 37 mm 	2	5	90	Bi-Color	36			9864-81
	2-0	3	150	Bi-Color	36			9829-51

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures


Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
1/2 CIRCLE 37 mm H 	0	3.5	75	Beige	36	•	9526-61	
1/2 CIRCLE 40 mm 	2-0	3	75	Bi-Color	36		9865-51	
	0	3.5	75	Bi-Color	36		9865-61	
	0	3.5	90	Bi-Color	36		9866-61	
	1	4	75	Bi-Color	36		9865-71	
	1	4	90	Bi-Color	36		9866-71	
	2	5	90	Bi-Color	36		9866-81	
1/2 CIRCLE 40 mm H 	0	3.5	75	Bi-Color	36		9868-61	
1/2 CIRCLE 48 mm 	2-0	3	75	Bi-Color	36		9867-51	
	0	3.5	75	Bi-Color	36		9867-61	
	1	4	75	Bi-Color	36		9867-71	
	1	4	90	Bi-Color	36		9885-71	
	2	5	75	Bi-Color	36		9867-81	
	2	5	90	Beige	36		9556-81	
	2	5	90	Bi-Color	36		9885-81	

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
1/2 CERCHIO 55 mm Super Rinforzato  HOS-14	2	5	75	Bi-Color	36		9850-81	
STRAIGHT 60 mm  SC-2	4-0	1.5	75	Beige	36		9108-31	
	3-0	2	75	Beige	36		9108-41	
	2-0	3	75	Beige	36		9108-51	

Cosmetic Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
3/8 CIRCLE 11 mm  P-10	6-0	0.7	45	Beige	36		9532-11	
	5-0	1	45	Beige	36		9532-21	
3/8 CIRCLE 13 mm  P-13	6-0	0.7	45	Beige	36		9536-11	
	5-0	1	45	Beige	36		9536-21	
	5-0	1	45	Beige	12		9536-23	
	4-0	1.5	45	Beige	36		9536-31	


Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Cosmetic Reverse Cutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
3/8 CIRCLE 16 mm 	5-0	1	45	Beige	36			9552-21
	4-0	1.5	75	Beige	36			9555-31
	3-0	2	75	Beige	36			9555-41
3/8 CIRCLE 19 mm 	5-0	1	45	Beige	36			9604-21
	5-0	1	75	Beige	36			9600-21
	4-0	1.5	45	Beige	36			9604-31
	4-0	1.5	45	Beige	12			9604-33
	4-0	1.5	75	Beige	36			9600-31
	3-0	2	45	Beige	36			9604-41
	3-0	2	45	Beige	12			9604-43
	3-0	2	75	Beige	36			9600-41
3/8 CIRCLE 24 mm 	4-0	1.5	45	Beige	36		•	9624-31
	4-0	1.5	75	Beige	36			9606-31
	3-0	2	45	Beige	36			9624-41
	3-0	2	75	Beige	36			9606-41
	2-0	3	45	Beige	36		•	9624-51
	2-0	3	75	Beige	36	SURGALLOY™		9606-51
1/2 CIRCLE 9 mm 	5-0	1	45	Beige	36			9704-21
1/2 CIRCLE 13 mm 	5-0	1	45	Beige	36			9712-21
	4-0	1.5	45	Beige	36			9712-31



Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Cosmetic Reverse Cutting							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
1/2 CIRCLE 16 mm 	4-0	1.5	45	Beige	36		9728-31



Spatula								
Description		Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
1/4 CIRCLE 8 mm		6-0	0.7	30	Beige	12	•	9560-13
		6-0	0.7	45	Beige	24		9509-12
		5-0	1	30	Beige	12	•	9560-23
		5-0	1	45	Beige	24		9509-22
Curvature	90°							
Chord Lgth	7.18 mm							
Radius	5.08 mm							

Taper Point							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
3/8 CIRCLE 17 mm 	4-0	1.5	75	Bi-Color	36		9812-31
	3-0	2	75	Bi-Color	36		9812-41
3/8 CIRCLE 22 mm 	4-0	1.5	75	Bi-Color	36		9583-31
	3-0	2	75	Bi-Color	36		9583-41




Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
3/8 CIRCLE 26 mm 	3-0	2	75	Bi-Color	36			9586-41
	2-0	3	75	Bi-Color	36			9586-51
3/8 CIRCLE 32 mm 	3-0	2	75	Bi-Color	36			9815-41
	2-0	3	75	Bi-Color	36			9815-51
	0	3.5	75	Bi-Color	36			9815-61
1/2 CIRCLE 13 mm 	6-0	0.7	75	Beige	36			9104-11
	5-0	1	75	Beige	36			9104-21
	4-0	1.5	75	Bi-Color	36			9810-31
	3-0	2	75	Bi-Color	36			9810-41
1/2 CIRCLE 17 mm 	6-0	0.7	75	Beige	36			9106-11
	5-0	1	75	Beige	36			9106-21
	4-0	1.5	5x45	Violet	36	D-TACH™	•	9109-37
	4-0	1.5	75	Beige	36			9106-31
	4-0	1.5	75	Bi-Color	36			9811-31
	3-0	2	75	Beige	36			9106-41
	3-0	2	75	Bi-Color	36			9811-41
	2-0	3	75	Bi-Color	36			9811-51
1/2 CIRCLE 20 mm 	5-0	1	75	Beige	36		•	9116-21
	4-0	1.5	5x45	Bi-Color	36	D-TACH™		9814-31
	4-0	1.5	75	Beige	36			9116-31
	4-0	1.5	75	Bi-Color	36			9813-31
	3-0	2	5x45	Bi-Color	36	D-TACH™		9814-41
	3-0	2	75	Beige	36			9116-41
	3-0	2	75	Bi-Color	36			9813-41
	3-0	2	75	Violet	36		•	9117-47
	2-0	3	75	Bi-Color	36			9813-51



Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 20 mm 	3-0	2	90	Bi-Color	36			9898-41
1/2 CIRCLE 22 mm 	4-0	1.5	75	Bi-Color	36			9227-31
	3-0	2	5x45	Bi-Color	36	D-TACH™		9789-41
	3-0	2	75	Bi-Color	36			9227-41
	3-0	2	5x75	Bi-Color	36	D-TACH™		9891-41
	2-0	3	5x45	Bi-Color	36	D-TACH™		9789-51
	2-0	3	75	Bi-Color	36			9227-51
	2-0	3	5x75	Bi-Color	36	D-TACH™		9891-51
1/2 CIRCLE 22 mm 	3-0	2	75	Bi-Color	36			9817-41
	2-0	3	75	Bi-Color	36			9817-51
1/2 CIRCLE 26 mm 	4-0	1.5	5x45	Beige	36	D-TACH™		9228-31
	4-0	1.5	5x45	Bi-Color	36	D-TACH™		9818-31
	4-0	1.5	75	Beige	36			9232-31
	4-0	1.5	75	Bi-Color	36			9819-31
	4-0	1.5	75	Violet	36		•	9233-37
	3-0	2	5x45	Beige	36	D-TACH™		9228-41
	3-0	2	5x45	Bi-Color	36	D-TACH™		9818-41
	3-0	2	75	Beige	36			9232-41
	3-0	2	75	Bi-Color	36			9819-41
	3-0	2	5x75	Bi-Color	36	D-TACH™		9917-41
	3-0	2	90	Bi-Color	36			9237-41
	2-0	3	5x45	Beige	36	D-TACH™		9228-51
	2-0	3	5x45	Bi-Color	36	D-TACH™		9818-51
	2-0	3	75	Beige	36			9232-51
	2-0	3	75	Bi-Color	36			9819-51
	2-0	3	5x75	Bi-Color	36	D-TACH™		9917-51
	2-0	3	90	Bi-Color	36			9237-51

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 26 mm 	0	3.5	75	Beige	36			9232-61
	0	3.5	75	Bi-Color	36			9819-61
	0	3.5	90	Bi-Color	36			9237-61
	1	4	75	Bi-Color	36			9819-71
1/2 CIRCLE 26 mm  STRAIGHT 44 mm 	3-0	2	90	Beige	36			9236-41
1/2 CIRCLE 26 mm 	3-0	2	90	Beige	36		•	9234-41
	3-0	2	90	Bi-Color	36			9820-41
	2-0	3	90	Bi-Color	36			9820-51
1/2 CIRCLE 27 mm 	3-0	2	75	Beige	36			9248-41
	3-0	2	75	Bi-Color	36			9822-41
	2-0	3	75	Beige	36			9248-51
	2-0	3	75	Bi-Color	36			9822-51
	0	3.5	5x45	Bi-Color	36	D-TACH™		9821-61
	0	3.5	75	Beige	36			9248-61
	0	3.5	75	Bi-Color	36			9822-61
	1	4	75	Bi-Color	36			9822-71
	1	4	5x75	Bi-Color	36	D-TACH™		9721-71


Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 30 mm 	3-0	2	75	Bi-Color	36			9825-41
	2-0	3	75	Beige	36			9855-51
	2-0	3	75	Bi-Color	36			9825-51
	0	3.5	75	Bi-Color	36			9825-61
	1	4	75	Bi-Color	36			9825-71
	1	4	90	Bi-Color	36			9924-71
1/2 CIRCLE 37 mm 	3-0	2	75	Bi-Color	36			9823-41
	2-0	3	75	Bi-Color	36			9823-51
	0	3.5	75	Bi-Color	36			9823-61
	1	4	75	Bi-Color	36			9823-71
1/2 CIRCLE 37 mm 	4-0	1.5	75	Beige	36			9266-31
	3-0	2	75	Beige	36			9266-41
	3-0	2	75	Bi-Color	36			9827-41
	3-0	2	90	Bi-Color	36			9828-41
	2-0	3	75	Beige	36			9266-51
	2-0	3	75	Bi-Color	36			9827-51
	2-0	3	75	Violet	36			9267-57
	2-0	3	90	Beige	36			9268-51
	2-0	3	90	Bi-Color	36			9828-51
	0	3.5	5x45	Bi-Color	36	D-TACH™		9830-61
	0	3.5	75	Beige	36			9266-61
	0	3.5	75	Bi-Color	36			9827-61
	0	3.5	90	Beige	36			9268-61
	0	3.5	90	Bi-Color	36			9828-61
	1	4	5x45	Bi-Color	36	D-TACH™		9830-71
	1	4	75	Beige	36			9266-71
	1	4	75	Bi-Color	36			9827-71
	1	4	90	Beige	36			9268-71
	1	4	90	Bi-Color	36			9828-71
	1	4	90	Violet	36			9269-77
	2	5	75	Bi-Color	36			9827-81
	2	5	90	Bi-Color	36			9828-81



Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 37 mm H  HGS-21	0	3.5	75	Bi-Color	36			9824-61
	0	3.5	90	Bi-Color	36			9841-61
	1	4	75	Bi-Color	36			9824-71
	1	4	90	Bi-Color	36			9841-71
	2	5	75	Beige	36			9258-81
	2	5	75	Bi-Color	36			9824-81
1/2 CIRCLE 40 mm  GS-24	3-0	2	75	Bi-Color	36			9832-41
	3-0	2	90	Bi-Color	36			9833-41
	2-0	3	75	Bi-Color	36			9832-51
	2-0	3	75	Violet	36		•	9285-57
	2-0	3	90	Bi-Color	36			9833-51
	0	3.5	75	Bi-Color	36			9832-61
	0	3.5	75	Violet	36		•	9285-67
	0	3.5	90	Bi-Color	36			9833-61
	1	4	5x45	Bi-Color	36	D-TACH™		9851-71
	1	4	5x60	Bi-Color	36	D-TACH™		9901-71
	1	4	75	Bi-Color	36			9832-71
	1	4	90	Beige	36			9286-71
	1	4	90	Bi-Color	36			9833-71
	2	5	5x45	Bi-Color	36	D-TACH™		9851-81
	2	5	5x60	Bi-Color	36	D-TACH™		9901-81
	2	5	75	Bi-Color	36			9832-81
	2	5	90	Bi-Color	36			9833-81




Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 48 mm  GS-25	2-0	3	75	Bi-Color	36			9835-51
	2-0	3	90	Bi-Color	36			9836-51
	0	3.5	75	Bi-Color	36			9835-61
	0	3.5	90	Beige	36			9304-61
	0	3.5	90	Bi-Color	36			9836-61
	1	4	5x45	Bi-Color	36	D-TACH™		9872-71
	1	4	5x60	Bi-Color	36	D-TACH™		9902-71
	1	4	75	Bi-Color	36			9835-71
	1	4	90	Bi-Color	36			9836-71
	2	5	5x45	Bi-Color	36	D-TACH™		9872-81
	2	5	5x60	Bi-Color	36	D-TACH™		9902-81
	2	5	75	Bi-Color	36			9835-81
	2	5	90	Bi-Color	36			9836-81
1/2 CIRCLE 60 mm  GS-34	2-0	3	75	Bi-Color	36			9837-51
	0	3.5	75	Bi-Color	36			9837-61
	1	4	75	Bi-Color	36			9837-71
1/2 CIRCLE 65 mm  GS-26	1	4	2x75	Beige	12			9314-73
	2	5	2x75	Beige	12			9314-83
	2	5	100	Beige	36			9322-81

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures


Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 76 mm 	3-0	2	75	Beige	36			9312-41
	2-0	3	75	Beige	36		•	9312-51
	0	3.5	75	Bi-Color	36		•	9840-61
	1	4	75	Bi-Color	36			9840-71
1/2 CIRCLE LOOPED 76 mm 	2	5	150	Bi-Color	24	150cm sut, 75cm looped		9848-82
5/8 CIRCLE 27 mm 	2-0	3	75	Bi-Color	36			9842-51
	0	3.5	75	Bi-Color	36			9842-61
5/8 CIRCLE 27 mm H 	0	3.5	45	Bi-Color	36			9849-61
	0	3.5	75	Bi-Color	36			9844-61
	1	4	75	Bi-Color	36			9844-71

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Taper Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
5/8 CIRCLE 37 mm 	2-0	3	75	Bi-Color	36		9843-51	
	0	3.5	75	Bi-Color	36		9843-61	
J NEEDLE 32 mm 	1	4	75	Bi-Color	12		9886-73	
STRAIGHT 44 mm 	2-0	3	75	Bi-Color	36	•	9845-51	

Blunt Taper								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
1/2 CIRCLE 37 mm 	2-0	3	90	Bi-Color	36		9042-51	
	0	3.5	90	Bi-Color	36		9042-61	

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures





Blunt Taper							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
1/2 CIRCLE 37 mm H  HBGS-21	2-0	3	90	Bi-Color	36		9035-51
	0	3.5	90	Bi-Color	36		9035-61
1/2 CIRCLE 40 mm  BGS-24	0	3.5	90	Bi-Color	36		9033-61
	1	4	90	Bi-Color	36		9033-71
1/2 CIRCLE 48 mm  BGS-25	0	3.5	90	Bi-Color	36		9031-61
	1	4	90	Bi-Color	36		9031-71
1/2 CIRCLE 65 mm  BGS-29	2	5	2x90	Bi-Color	36	PROTECT-POINT™	9037-81

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Blunt Point								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
3/8 CIRCLE 64 mm 	0	3.5	90	Beige	36	PROTECT-POINT™		9142-61
	1	4	90	Beige	36	PROTECT-POINT™		9142-71

Tapercutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 22 mm 	1	4	75	Beige	36			9252-71
1/2 CIRCLE 26 mm 	3-0	2	75	Bi-Color	36			9251-41
	2-0	3	75	Bi-Color	36			9251-51
1/2 CIRCLE 27 mm 	2-0	3	90	Beige	36		•	9253-51
	0	3.5	90	Beige	36			9253-61


Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Tapercutting								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 37 mm  KV-34	3-0	2	90	Bi-Color	36			9831-41
	2-0	3	90	Beige	36			9857-51
	2-0	3	90	Bi-Color	36			9831-51
	0	3.5	90	Bi-Color	36			9831-61
	1	4	90	Beige	36			9857-71
	1	4	90	Bi-Color	36			9831-71
	2	5	90	Bi-Color	36			9831-81
1/2 CIRCLE 40 mm  KV-37	0	3.5	75	Bi-Color	36			9834-61
	2	5	90	Bi-Color	36			9930-81
1/2 CIRCLE 48 mm  KV-40	0	3.5	90	Bi-Color	36			9847-61
	1	4	90	Bi-Color	36			9847-71
	2	5	90	Bi-Color	36			9847-81
1/2 CIRCLE LOOPED 48 mm  KV-40	1	4	150	Bi-Color	24	150cm sut, 75cm looped		9846-72
	2	5	150	Bi-Color	24	150cm sut, 75cm looped		9846-82



Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Multipurpose								
Description		Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
3/8 CIRCLE 13 mm		5-0	1	45	Beige	36	•	6313-84
3/8 CIRCLE 16 mm		5-0	1	45	Beige	36	•	6316-84
		4-0	1.5	45	Beige	36		6316-85
		3-0	2	45	Beige	36	•	6316-86
3/8 CIRCLE 35 mm		2-0	3	75	Beige	36	•	6335-87
1/2 CIRCLE 16 mm		4-0	1.5	75	Beige	36	•	6416-85
1/2 CIRCLE 20 mm		4-0	1.5	75	Bi-Color	36		6420-75
		3-0	2	75	Beige	36	•	6420-86
		3-0	2	75	Bi-Color	36	•	6420-76




Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Multipurpose								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 20 mm 	2-0	3	75	Bi-Color	36			6421-77
	0	3.5	75	Bi-Color	36		•	6421-78
1/2 CIRCLE 25 mm 	3-0	2	75	Beige	36		•	6425-86
	3-0	2	75	Bi-Color	36		•	6425-76
	2-0	3	75	Bi-Color	36			6425-77
1/2 CIRCLE 30 mm 	3-0	2	75	Bi-Color	36		•	6430-76
	2-0	3	75	Bi-Color	36			6430-77
1/2 CIRCLE 30 mm 	2-0	3	75	Bi-Color	36			6431-77
	0	3.5	75	Bi-Color	36			6431-78
1/2 CIRCLE 35 mm 	1	4	75	Bi-Color	36		•	6435-79
1/2 CIRCLE 40 mm 	2	5	90	Bi-Color	36			6440-70

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Multipurpose								
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code	
1/2 CIRCLE 40 mm 	2-0	3	75	Beige	36	•	6441-87	
	2-0	3	75	Bi-Color	36		6441-77	
	0	3.5	75	Beige	36	•	6441-88	
	0	3.5	75	Bi-Color	36		6441-78	
	1	4	75	Beige	36		6441-89	
	1	4	75	Bi-Color	36		6441-79	
1/2 CIRCLE 45 mm 	1	4	75	Bi-Color	36		6445-79	
1/2 CIRCLE 50 mm 	2	5	90	Bi-Color	36		6450-70	
1/2 CIRCLE 50 mm 	2-0	3	75	Bi-Color	36		6451-77	
	0	3.5	75	Bi-Color	36		6451-78	
	1	4	75	Bi-Color	36		6451-79	

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Multipurpose							
Description	Size	Size (Metric)	Length (cm)	Color	Units per box	Phasing-out	Reorder code
5/8 CIRCLE 110 mm  5.99 	2	5	150	Bi-Color	12	•	6599-70
STRAIGHT 45 mm  1.45 	4-0	1.5	75	Beige	36	•	6145-85
	3-0	2	75	Beige	36	•	6145-86
STRAIGHT 75 mm  1.75 	2-0	3	75	Beige	12		6175-87

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Index					
Reorder code	Page	Reorder code	Page	Reorder code	Page
6145-85	24	9008-42	2	9228-51	11
6145-86	24	9008-52	2	9232-31	11
6175-87	24	9008-62	2	9232-41	11
6313-84	21	9017-32	1	9232-51	11
6316-84	21	9017-42	1	9232-61	12
6316-85	21	9017-52	1	9233-37	11
6316-86	21	9017-62	1	9234-41	12
6335-87	21	9017-72	1	9236-41	12
6416-85	21	9017-82	1	9237-41	11
6420-75	21	9031-61	18	9237-51	11
6420-76	21	9031-71	18	9237-61	12
6420-86	21	9033-61	18	9248-41	12
6421-77	22	9033-71	18	9248-51	12
6421-78	22	9035-51	18	9248-61	12
6425-76	22	9035-61	18	9251-41	19
6425-77	22	9037-81	18	9251-51	19
6425-86	22	9042-51	17	9252-71	19
6430-76	22	9042-61	17	9253-51	19
6430-77	22	9104-11	10	9253-61	19
6431-77	22	9104-21	10	9258-81	14
6431-78	22	9106-11	10	9266-31	13
6435-79	22	9106-21	10	9266-41	13
6440-70	22	9106-31	10	9266-51	13
6441-77	23	9106-41	10	9266-61	13
6441-78	23	9108-31	7	9266-71	13
6441-79	23	9108-41	7	9267-57	13
6441-87	23	9108-51	7	9268-51	13
6441-88	23	9109-37	10	9268-61	13
6441-89	23	9116-21	10	9268-71	13
6445-79	23	9116-31	10	9269-77	13
6450-70	23	9116-41	10	9285-57	14
6451-77	23	9117-47	10	9285-67	14
6451-78	23	9142-61	19	9286-71	14
6451-79	23	9142-71	19	9304-61	15
6599-70	24	9227-31	11	9312-41	16
9004-42	1	9227-41	11	9312-51	16
9004-52	1	9227-51	11	9314-73	15
9004-62	1	9228-31	11	9314-83	15
9004-72	1	9228-41	11	9322-81	15

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Index					
Reorder code	Page	Reorder code	Page	Reorder code	Page
9509-12	9	9606-41	8	9801-61	2
9509-22	9	9606-51	8	9801-71	2
9520-51	5	9608-21	3	9801-81	2
9526-61	6	9608-31	3	9802-31	1
9532-11	7	9608-41	3	9802-41	1
9532-21	7	9612-21	3	9802-51	1
9536-11	7	9612-31	3	9802-61	1
9536-21	7	9612-41	3	9802-71	1
9536-23	7	9618-31	3	9803-31	1
9536-31	7	9618-41	3	9803-41	1
9538-11	3	9618-51	3	9803-51	1
9538-21	3	9619-31	3	9803-61	1
9538-31	3	9619-41	3	9803-71	1
9552-21	8	9619-51	3	9803-81	1
9555-31	8	9624-31	8	9804-33	2
9555-41	8	9624-41	8	9804-43	2
9556-81	6	9624-51	8	9804-53	2
9560-13	9	9626-41	4	9804-63	2
9560-23	9	9626-51	4	9804-73	2
9583-31	9	9626-61	4	9805-31	1
9583-41	9	9626-71	4	9805-41	1
9584-41	5	9627-51	4	9805-51	1
9584-51	5	9627-61	4	9805-61	1
9584-61	5	9628-51	4	9805-73	1
9584-71	5	9641-57	3	9805-83	1
9586-41	10	9704-21	8	9806-41	1
9586-51	10	9712-21	8	9806-51	1
9600-21	8	9712-31	8	9807-43	2
9600-31	8	9721-71	12	9807-53	2
9600-41	8	9728-31	9	9807-63	2
9602-21	3	9735-41	4	9808-33	2
9602-31	3	9740-41	4	9808-37	2
9602-41	3	9740-51	4	9808-43	2
9604-21	8	9746-51	4	9808-47	2
9604-31	8	9789-41	11	9808-53	2
9604-33	8	9789-51	11	9808-57	2
9604-41	8	9801-31	2	9808-63	2
9604-43	8	9801-41	2	9808-67	2
9606-31	8	9801-51	2	9808-73	2

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Index					
Reorder code	Page	Reorder code	Page	Reorder code	Page
9808-77	2	9822-71	12	9833-81	14
9808-83	2	9823-41	13	9834-61	20
9809-32	1	9823-51	13	9835-51	15
9809-42	1	9823-61	13	9835-61	15
9809-52	1	9823-71	13	9835-71	15
9809-62	1	9824-61	14	9835-81	15
9809-72	1	9824-71	14	9836-51	15
9809-82	1	9824-81	14	9836-61	15
9810-31	10	9825-41	13	9836-71	15
9810-41	10	9825-51	13	9836-81	15
9811-31	10	9825-61	13	9837-51	15
9811-41	10	9825-71	13	9837-61	15
9811-51	10	9827-41	13	9837-71	15
9812-31	9	9827-51	13	9839-51	3
9812-41	9	9827-61	13	9840-61	16
9813-31	10	9827-71	13	9840-71	16
9813-41	10	9827-81	13	9841-61	14
9813-51	10	9828-41	13	9841-71	14
9814-31	10	9828-51	13	9842-51	16
9814-41	10	9828-61	13	9842-61	16
9815-41	10	9828-71	13	9843-51	17
9815-51	10	9828-81	13	9843-61	17
9815-61	10	9829-51	5	9844-61	16
9817-41	11	9830-61	13	9844-71	16
9817-51	11	9830-71	13	9845-51	17
9818-31	11	9831-41	20	9846-72	20
9818-41	11	9831-51	20	9846-82	20
9818-51	11	9831-61	20	9847-61	20
9819-31	11	9831-71	20	9847-71	20
9819-41	11	9831-81	20	9847-81	20
9819-51	11	9832-41	14	9848-82	16
9819-61	12	9832-51	14	9849-61	16
9819-71	12	9832-61	14	9850-81	7
9820-41	12	9832-71	14	9851-71	14
9820-51	12	9832-81	14	9851-81	14
9821-61	12	9833-41	14	9855-51	13
9822-41	12	9833-51	14	9857-51	20
9822-51	12	9833-61	14	9857-71	20
9822-61	12	9833-71	14	9860-31	4

Dexon™ II Polyglycolic Acid, Coated, Braided Synthetic Absorbable Sutures

Index					
Reorder code	Page	Reorder code	Page	Reorder code	Page
9861-41	4	9898-41	11		
9861-51	4	9901-71	14		
9861-61	4	9901-81	14		
9862-61	5	9902-71	15		
9862-71	5	9902-81	15		
9863-41	5	9917-41	11		
9863-51	5	9917-51	11		
9863-61	5	9924-71	13		
9863-71	5	9930-81	20		
9863-81	5				
9864-51	5				
9864-61	5				
9864-71	5				
9864-81	5				
9865-51	6				
9865-61	6				
9865-71	6				
9866-61	6				
9866-71	6				
9866-81	6				
9867-51	6				
9867-61	6				
9867-71	6				
9867-81	6				
9868-61	6				
9872-71	15				
9872-81	15				
9875-41	4				
9875-51	4				
9877-42	1				
9877-52	1				
9877-62	1				
9877-72	1				
9877-82	2				
9885-71	6				
9885-81	6				
9886-73	17				
9891-41	11				
9891-51	11				