

Interventional
imaging

▶ OCTOFLUX

Procedures made easier with sets
designed for injectors

- Filling and injection sets
- For single head and double head injectors
- Patient lines with non-return valves



PEROUSE
MEDICAL



FILLING AND INJECTION SET (maximum pressure 300 PSI)

	References	Description	Packaging
	0170NF	Single head injector: <ul style="list-style-type: none"> • 2 non-return valves - contrast media line 100 cm with a drip chamber - patient line 100 cm 	25 unit-box
	0170PS	Double head injector: <ul style="list-style-type: none"> • 4 non-return valves - 1 contrast media line 100 cm with a drip chamber + tubing of 10 cm - 1 saline line 100 cm with a drip chamber + tubing of 40 cm - 1 patient line 10 cm 	25 unit-box
	0170NQ	Double head injector: <ul style="list-style-type: none"> - patient line 150 cm - syringe to syringe line 50 cm 	25 unit-box
	0170RK	Double head injector: <ul style="list-style-type: none"> • 1 non-return valve - patient line 150 cm - syringe to syringe line 25 cm 	50 unit-box
	0170SR	Double head injector: <ul style="list-style-type: none"> • 2 non-return valves - patient line 150 cm - syringe to syringe line 25 cm 	50 unit-box

PATIENT LINE (maximum pressure 300PSI)

	References	Description	Packaging
	0170NB	<ul style="list-style-type: none"> • 1 non-return valve - patient line 20 cm 	50 unit-box
	0170NG	<ul style="list-style-type: none"> • 1 non-return valve - patient line 40 cm 	50 unit-box
	0170NN	<ul style="list-style-type: none"> • 1 non-return valve - patient line 60 cm 	50 unit-box
	0170NT	<ul style="list-style-type: none"> • 1 non-return valve - patient line 120 cm 	50 unit-box

