

VASCULAR ACCESS  
Accessories for implantable ports



Huber  
**needles**



Value Life

# Table of Content

## **Non-safety huber needles with tubing**

Hubsite 4

Polyperf 6

## **Non-safety huber needles without tubing**

Huber needle 8

Inject-site 9

## **Safety huber needles with manual positive pressure**

Perfusafe 2 11

PPS® 14

## **Safety huber needles with automatic positive pressure**

PPS® Flow+ 17



Non-safety

## Non-safety huber needles with tubing

- Hubsite



### Hubsite

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
I239.715	22 G	0.7 mm	15 mm	yes	10 / carton box
I239.720	22 G	0.7 mm	20 mm	yes	10 / carton box
I239.725	22 G	0.7 mm	25 mm	yes	10 / carton box
I239.730	22 G	0.7 mm	30 mm	yes	10 / carton box
I239.915	20 G	0.9 mm	15 mm	yes	10 / carton box
I239.920	20 G	0.9 mm	20 mm	yes	10 / carton box
I239.925	20 G	0.9 mm	25 mm	yes	10 / carton box
I239.930	20 G	0.9 mm	30 mm	yes	10 / carton box
I239.115	19 G	0.7 mm	15 mm	yes	10 / carton box
I239.120	19 G	0.7 mm	20 mm	yes	10 / carton box
I239.125	19 G	0.7 mm	25 mm	yes	10 / carton box
I239.130	19 G	0.7 mm	30 mm	yes	10 / carton box

### Hubsite with lateral port and autoflush

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
I251.715	22 G	0.7 mm	15 mm	yes	10 / carton box
I251.720	22 G	0.7 mm	20 mm	yes	10 / carton box
I251.725	22 G	0.7 mm	25 mm	yes	10 / carton box
I251.730	22 G	0.7 mm	30 mm	yes	10 / carton box
I251.915	20 G	0.9 mm	15 mm	yes	10 / carton box
I251.920	20 G	0.9 mm	20 mm	yes	10 / carton box
I251.925	20 G	0.9 mm	25 mm	yes	10 / carton box
I251.930	20 G	0.9 mm	30 mm	yes	10 / carton box
I251.115	19 G	1.1 mm	15 mm	yes	10 / carton box
I251.120	19 G	1.1 mm	20 mm	yes	10 / carton box
I251.125	19 G	1.1 mm	25 mm	yes	10 / carton box
I251.130	19 G	1.1 mm	30 mm	yes	10 / carton box

- Polyperf



### Polyperf without lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE581507*	22 G	1.1 mm	15 mm	no	10 / carton box
VPE581707	22 G	1.1 mm	17 mm	no	10 / carton box
VPE582007	22 G	1.1 mm	20 mm	no	10 / carton box
VPE582507	22 G	1.1 mm	25 mm	no	10 / carton box
VPE583007*	22 G	1.1 mm	30 mm	no	10 / carton box
VPE583507*	22 G	1.1 mm	35 mm	no	10 / carton box
VPE581509*	20 G	0.9 mm	15 mm	no	10 / carton box
VPE581709	20 G	0.9 mm	17 mm	no	10 / carton box
VPE582009	20 G	0.9 mm	20 mm	no	10 / carton box
VPE582509	20 G	0.9 mm	25 mm	no	10 / carton box
VPE583009*	20 G	0.9 mm	30 mm	no	10 / carton box
VPE583509*	20 G	0.9 mm	35 mm	no	10 / carton box
VPE581511	19 G	1.1 mm	15 mm	no	10 / carton box
VPE581711	19 G	1.1 mm	17 mm	no	10 / carton box
VPE582011	19 G	1.1 mm	20 mm	no	10 / carton box
VPE582511	19 G	1.1 mm	25 mm	no	10 / carton box
VPE583011*	19 G	1.1 mm	30 mm	no	10 / carton box
VPE583511*	19 G	1.1 mm	35 mm	no	10 / carton box

\*These codes are only available upon request (minimum lead time is 3 months).

### Polyperf with lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE592007	22 G	0.7 mm	20 mm	no	10 / carton box
VPE592507*	22 G	0.7 mm	25 mm	no	10 / carton box
VPE592009	20 G	0.9 mm	20 mm	no	10 / carton box
VPE592509	20 G	0.9 mm	25 mm	no	10 / carton box
VPE592011	19 G	1.1 mm	20 mm	no	10 / carton box
VPE592511	19 G	1.1 mm	25 mm	no	10 / carton box

\*These codes are only available upon request (minimum lead time is 3 months).

## Non-safety huber needles without tubing

- Huber needle



### Huber needle with straight bevel

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE512506	23 G	0.6 mm	25 mm	no	50 / carton box
VPE512511*	22 G	0.7 mm	25 mm	no	50 / carton box
VPE513811*	22 G	0.7 mm	38 mm	no	50 / carton box
VPE512509	20 G	0.9 mm	25 mm	no	50 / carton box
VPE512507	19 G	1.1 mm	25 mm	no	50 / carton box
VPE513807	19 G	1.1 mm	38 mm	no	50 / carton box

\*These codes are only available upon request (minimum lead time is 3 months).

### Huber needle with angled bevel

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE522506*	23 G	0.6 mm	25 mm	no	50 / carton box
VPE522507	22 G	0.7 mm	25 mm	no	50 / carton box
VPE523807*	22 G	0.7 mm	38 mm	no	50 / carton box
VPE522509	20 G	0.9 mm	25 mm	no	50 / carton box
VPE523809	20 G	0.9 mm	38 mm	no	50 / carton box
VPE522511	19 G	1.1 mm	25 mm	no	50 / carton box
VPE523811*	19 G	1.1 mm	38 mm	no	50 / carton box

\*These codes are only available upon request (minimum lead time is 3 months).

- Inject-site



### Inject-site with straight bevel

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
1057.27	22 G	0.7 mm	25 mm	no	30 / carton box
1057.29	20 G	0.9 mm	25 mm	no	30 / carton box

### Inject-site with angled bevel

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
1058.07	22 G	0.7 mm	20 mm	no	30 / carton box
1058.27	22 G	0.7 mm	25 mm	no	30 / carton box
1058.09	20 G	0.9 mm	20 mm	no	30 / carton box
1058.29	20 G	0.9 mm	25 mm	no	30 / carton box


# Safety

## Safety huber needles with manual positive pressure

- Perfusafe 2

Perfusafe2 is a safety infusion device with a Huber bevel for use with implantable ports.

Perfusafe2 prevents needlestick injuries caused by the needle rebound effect. The needle is automatically directed into a safety cage offering protection from accidental needlestick injury, with no change in the handling technique.



### Perfusafe 2 without lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
5249.715	22 G	0.7 mm	15 mm	yes	10 / carton box
5249.720	22 G	0.7 mm	20 mm	yes	10 / carton box
5249.725	22 G	0.7 mm	25 mm	yes	10 / carton box
5249.730	22 G	0.7 mm	30 mm	yes	10 / carton box
5249.915	20 G	0.9 mm	15 mm	yes	10 / carton box
5249.920	20 G	0.9 mm	20 mm	yes	10 / carton box
5249.925	20 G	0.9 mm	25 mm	yes	10 / carton box
5249.930	20 G	0.9 mm	30 mm	yes	10 / carton box
5249.115	19 G	1.1 mm	15 mm	yes	10 / carton box
5249.120	19 G	1.1 mm	20 mm	yes	10 / carton box
5249.125	19 G	1.1 mm	25 mm	yes	10 / carton box
5249.130	19 G	1.1 mm	30 mm	yes	10 / carton box

### Perfusafe 2 (without lateral port) with autoflush

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
5249.916	22 G	0.7 mm	15 mm	yes	10 / carton box
5249.921	22 G	0.7 mm	20 mm	yes	10 / carton box
5249.926	22 G	0.7 mm	25 mm	yes	10 / carton box
5249.916	20 G	0.9 mm	15 mm	yes	10 / carton box
5249.921	20 G	0.9 mm	20 mm	yes	10 / carton box
5249.926	20 G	0.9 mm	25 mm	yes	10 / carton box
5249.116	19 G	1.1 mm	15 mm	yes	10 / carton box
5249.121	19 G	1.1 mm	20 mm	yes	10 / carton box
5249.126	19 G	1.1 mm	25 mm	yes	10 / carton box

### Perfusafe 2 (without lateral port) with bionector

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
5249.717	22 G	0.7 mm	15 mm	yes	10 / carton box
5249.722	22 G	0.7 mm	20 mm	yes	10 / carton box
5249.727	22 G	0.7 mm	25 mm	yes	10 / carton box
5249.732	22 G	0.7 mm	30 mm	yes	10 / carton box
5249.917	20 G	0.9 mm	15 mm	yes	10 / carton box
5249.922	20 G	0.9 mm	20 mm	yes	10 / carton box
5249.927	20 G	0.9 mm	25 mm	yes	10 / carton box
5249.932	20 G	0.9 mm	30 mm	yes	10 / carton box
5249.117	19 G	1.1 mm	15 mm	yes	10 / carton box
5249.122	19 G	1.1 mm	20 mm	yes	10 / carton box
5249.127	19 G	1.1 mm	25 mm	yes	10 / carton box
5249.132	19 G	1.1 mm	30 mm	yes	10 / carton box

### Perfusafe 2 with lateral port and autoflush

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
5250.715	22 G	0.7 mm	15 mm	yes	10 / carton box
5250.720	22 G	0.7 mm	20 mm	yes	10 / carton box
5250.725	22 G	0.7 mm	25 mm	yes	10 / carton box
5250.730	22 G	0.7 mm	30 mm	yes	10 / carton box
5250.915	20 G	0.9 mm	15 mm	yes	10 / carton box
5250.920	20 G	0.9 mm	20 mm	yes	10 / carton box
5250.925	20 G	0.9 mm	25 mm	yes	10 / carton box
5250.930	20 G	0.9 mm	30 mm	yes	10 / carton box
5250.115	19 G	1.1 mm	15 mm	yes	10 / carton box
5250.120	19 G	1.1 mm	20 mm	yes	10 / carton box
5250.125	19 G	1.1 mm	25 mm	yes	10 / carton box
5250.130	19 G	1.1 mm	30 mm	yes	10 / carton box

### Perfusafe 2 with lateral port and bionector

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
5250.717	22 G	0.7 mm	15 mm	yes	10 / carton box
5250.722	22 G	0.7 mm	20 mm	yes	10 / carton box
5250.727	22 G	0.7 mm	25 mm	yes	10 / carton box
5250.732	22 G	0.7 mm	30 mm	yes	10 / carton box
5250.917	20 G	0.9 mm	15 mm	yes	10 / carton box
5250.922	20 G	0.9 mm	20 mm	yes	10 / carton box
5250.927	20 G	0.9 mm	25 mm	yes	10 / carton box
5250.932	20 G	0.9 mm	30 mm	yes	10 / carton box
5250.117	19 G	1.1 mm	15 mm	yes	10 / carton box
5250.122	19 G	1.1 mm	20 mm	yes	10 / carton box
5250.127	19 G	1.1 mm	25 mm	yes	10 / carton box
5250.132	19 G	1.1 mm	30 mm	yes	10 / carton box

### QimoPerfusafe

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
7242.720	22 G	0.7 mm	20 mm	yes	10 / carton box
7242.725	22 G	0.7 mm	25 mm	yes	10 / carton box
7242.920	20 G	0.9 mm	20 mm	yes	10 / carton box
7242.925	20 G	0.9 mm	25 mm	yes	10 / carton box

• PPS®

PPS are safety infusion device with a Huber bevel for use with implantable ports. PPS prevents needlestick injuries caused by the needle rebound effect. The Withdrawal of the needle is possible with only one hand.

PPS® without lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE601507	22 G	0.7 mm	15 mm	no	12 / carton box
VPE601707	22 G	0.7 mm	17 mm	no	12 / carton box
VPE602007	22 G	0.7 mm	20 mm	no	12 / carton box
VPE602507	22 G	0.7 mm	25 mm	no	12 / carton box
VPE603007*	22 G	0.7 mm	30 mm	no	12 / carton box
VPE603507*	22 G	0.7 mm	35 mm	no	12 / carton box
VPE601509	20 G	0.9 mm	15 mm	no	12 / carton box
VPE601709	20 G	0.9 mm	17 mm	no	12 / carton box
VPE602009	20 G	0.9 mm	20 mm	no	12 / carton box
VPE602509	20 G	0.9 mm	25 mm	no	12 / carton box
VPE603009*	20 G	0.9 mm	30 mm	no	12 / carton box
VPE603509*	20 G	0.9 mm	35 mm	no	12 / carton box
VPE601511	19 G	1.1 mm	15 mm	no	12 / carton box
VPE601711	19 G	1.1 mm	17 mm	no	12 / carton box
VPE602011	19 G	1.1 mm	20 mm	no	12 / carton box
VPE602511	19 G	1.1 mm	25 mm	no	12 / carton box
VPE603011*	19 G	1.1 mm	30 mm	no	12 / carton box
VPE603511*	19 G	1.1 mm	35 mm	no	12 / carton box

PPS® with lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE611507	22 G	0.7 mm	15 mm	no	12 / carton box
VPE611707	22 G	0.7 mm	17 mm	no	12 / carton box
VPE612007	22 G	0.7 mm	20 mm	no	12 / carton box
VPE612507	22 G	0.7 mm	25 mm	no	12 / carton box
VPE613007	22 G	0.7 mm	30 mm	no	12 / carton box
VPE613507	22 G	0.7 mm	35 mm	no	12 / carton box
VPE611509	20 G	0.9 mm	15 mm	no	12 / carton box
VPE611709	20 G	0.9 mm	17 mm	no	12 / carton box
VPE612009	20 G	0.9 mm	20 mm	no	12 / carton box
VPE612509	20 G	0.9 mm	25 mm	no	12 / carton box
VPE613009	20 G	0.9 mm	30 mm	no	12 / carton box
VPE613509	20 G	0.9 mm	35 mm	no	12 / carton box

PPS® with lateral port (continuation)

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE611511	19 G	1.1 mm	15 mm	no	12 / carton box
VPE611711	19 G	1.1 mm	17 mm	no	12 / carton box
VPE612011	19 G	1.1 mm	20 mm	no	12 / carton box
VPE612511	19 G	1.1 mm	25 mm	no	12 / carton box
VPE613011	19 G	1.1 mm	30 mm	no	12 / carton box
VPE613511	19 G	1.1 mm	35 mm	no	12 / carton box

PPS® CT Rated without lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE801507	22 G	0.7 mm	15 mm	yes	12 / carton box
VPE801707	22 G	0.7 mm	17 mm	yes	12 / carton box
VPE802007	22 G	0.7 mm	20 mm	yes	12 / carton box
VPE802507	22 G	0.7 mm	25 mm	yes	12 / carton box
VPE803007	22 G	0.7 mm	30 mm	yes	12 / carton box
VPE803507	22 G	0.7 mm	35 mm	yes	12 / carton box
VPE801509	20 G	0.9 mm	15 mm	yes	12 / carton box
VPE801709	20 G	0.9 mm	17 mm	yes	12 / carton box
VPE802009	20 G	0.9 mm	20 mm	yes	12 / carton box
VPE802509	20 G	0.9 mm	25 mm	yes	12 / carton box
VPE803009	20 G	0.9 mm	30 mm	yes	12 / carton box
VPE803509	20 G	0.9 mm	35 mm	yes	12 / carton box
VPE801511	19 G	1.1 mm	15 mm	yes	12 / carton box
VPE801711	19 G	1.1 mm	17 mm	yes	12 / carton box
VPE802011	19 G	1.1 mm	20 mm	yes	12 / carton box
VPE802511	19 G	1.1 mm	25 mm	yes	12 / carton box
VPE803011	19 G	1.1 mm	30 mm	yes	12 / carton box
VPE803511	19 G	1.1 mm	35 mm	yes	12 / carton box

PPS® CT Rated with lateral port




Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE811507	22 G	0.7 mm	15 mm	yes	12 / carton box
VPE811707	22 G	0.7 mm	17 mm	yes	12 / carton box
VPE812007	22 G	0.7 mm	20 mm	yes	12 / carton box
VPE812507	22 G	0.7 mm	25 mm	yes	12 / carton box
VPE813007	22 G	0.7 mm	30 mm	yes	12 / carton box
VPE813507	22 G	0.7 mm	35 mm	yes	12 / carton box
VPE811509	20 G	0.9 mm	15 mm	yes	12 / carton box
VPE811709	20 G	0.9 mm	17 mm	yes	12 / carton box
VPE812009	20 G	0.9 mm	20 mm	yes	12 / carton box
VPE812509	20 G	0.9 mm	25 mm	yes	12 / carton box
VPE813009	20 G	0.9 mm	30 mm	yes	12 / carton box
VPE813509	20 G	0.9 mm	35 mm	yes	12 / carton box
VPE811511	19 G	1.1 mm	15 mm	yes	12 / carton box
VPE811711	19 G	1.1 mm	17 mm	yes	12 / carton box
VPE812011	19 G	1.1 mm	20 mm	yes	12 / carton box
VPE812511	19 G	1.1 mm	25 mm	yes	12 / carton box
VPE813011	19 G	1.1 mm	30 mm	yes	12 / carton box
VPE813511	19 G	1.1 mm	35 mm	yes	12 / carton box

Safety huber needles with automatic positive pressure

- PPS® Flow+

PPS Flow+ is a new generation of Safety Huber needles with automatic positive pressure .

PPS® Flow+ without lateral port



Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE701507	22 G	0.7 mm	15 mm	no	12 / carton box
VPE701707	22 G	0.7 mm	17 mm	no	12 / carton box
VPE702007	22 G	0.7 mm	20 mm	no	12 / carton box
VPE702507	22 G	0.7 mm	25 mm	no	12 / carton box
VPE703007	22 G	0.7 mm	30 mm	no	12 / carton box
VPE703507	22 G	0.7 mm	35 mm	no	12 / carton box
VPE701509	20 G	0.9 mm	15 mm	no	12 / carton box
VPE701709	20 G	0.9 mm	17 mm	no	12 / carton box
VPE702009	20 G	0.9 mm	20 mm	no	12 / carton box
VPE702509	20 G	0.9 mm	25 mm	no	12 / carton box
VPE703009	20 G	0.9 mm	30 mm	no	12 / carton box
VPE703509	20 G	0.9 mm	35 mm	no	12 / carton box
VPE701511	19 G	1.1 mm	15 mm	no	12 / carton box
VPE701711	19 G	1.1 mm	17 mm	no	12 / carton box
VPE702011	19 G	1.1 mm	20 mm	no	12 / carton box
VPE702511	19 G	1.1 mm	25 mm	no	12 / carton box
VPE703011	19 G	1.1 mm	30 mm	no	12 / carton box
VPE703511	19 G	1.1 mm	35 mm	no	12 / carton box




PPS® Flow+ with lateral port

Code	Gauge	Ext. Ø	Length	CT Rated	Packaging
VPE711507	22 G	0.7 mm	15 mm	no	12 / carton box
VPE711707	22 G	0.7 mm	17 mm	no	12 / carton box
VPE712007	22 G	0.7 mm	20 mm	no	12 / carton box
VPE712507	22 G	0.7 mm	25 mm	no	12 / carton box
VPE713007	22 G	0.7 mm	30 mm	no	12 / carton box
VPE713507	22 G	0.7 mm	35 mm	no	12 / carton box
VPE711509	20 G	0.9 mm	15 mm	no	12 / carton box
VPE711709	20 G	0.9 mm	17 mm	no	12 / carton box
VPE712009	20 G	0.9 mm	20 mm	no	12 / carton box
VPE712509	20 G	0.9 mm	25 mm	no	12 / carton box
VPE713009	20 G	0.9 mm	30 mm	no	12 / carton box
VPE713509	20 G	0.9 mm	35 mm	no	12 / carton box
VPE711511	19 G	1.1 mm	15 mm	no	12 / carton box
VPE711711	19 G	1.1 mm	17 mm	no	12 / carton box
VPE712011	19 G	1.1 mm	20 mm	no	12 / carton box
VPE712511	19 G	1.1 mm	25 mm	no	12 / carton box
VPE713011	19 G	1.1 mm	30 mm	no	12 / carton box
VPE713511	19 G	1.1 mm	35 mm	no	12 / carton box

For further information, please contact: [marketingbenelux@vygon.com](mailto:marketingbenelux@vygon.com)

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 Chaussée de Haecht • 1130 BRUSSELS • BELGIUM  
 Tel: +32 (0)2 706 09 50 • [www.vygon.be](http://www.vygon.be)  
 Vygon Nederland B.V. • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • NEDERLAND  
 Tel: +31 (0)40 208 93 80 • [www.vygon.nl](http://www.vygon.nl)  
 VY16076 •  Vygon Benelux

