

Lames de laryngoscope

à usage unique




design
amélioré

Lames de laryngoscope à usage unique

Introduction

L'efficacité des procédures de nettoyage et de désinfection des lames de laryngoscope préoccupe aujourd'hui les acteurs de la santé. De nombreux articles cliniques apportent des éléments qui confortent ces interrogations^{1, 2, 3, 4, 5, 6}. La lame à usage unique libère des contraintes et des limites de la désinfection et de la stérilisation, ce qui constitue une avancée réelle du point de vue de la sécurité sanitaire pour les malades et les utilisateurs.

Nos lames, Miller tailles 0 et 1, et MacIntosh tailles 2, 3 et 4, identifiées par le logo vert Intersurgical sur le talon, sont compatibles avec tous les manches à bague verte ISO7376-3. Elles apportent une solution économique au problème des infections croisées.



Lame de laryngoscope Miller - taille 0

Code: 7040

Boîte de: 8



Lame de laryngoscope Miller - taille 1

Code: 7041

Boîte de: 8



Lame de laryngoscope Macintosh - taille 2

Code: 7042

Boîte de: 24



Lame de laryngoscope Macintosh - taille 3

Code: 7043

Boîte de: 20



Lame de laryngoscope Macintosh - taille 4

Code: 7044

Boîte de: 20

Références

1. **'Decontamination of Laryngoscopes: is current practice effective?'**
Andrew Hutchinson – Princess Royal Hospital, Telford (Technic Issue 223 March 2002)
2. **'Decontamination of Laryngoscopes: a survey of national practice'**
M.D.Esler et al – Royal London Hospital, London (Anaesthesia 1999, 54, pages 587-592)
3. **'Incidence of visible and occult blood on Lame de laryngoscopes and handles'**
R. Phillips & W.P. Monaghan – Fort Worth, Texas (Journal of the American Association of Nurse Anaesthetists June 1997, Vol 65, No 3, page 241-246)
4. **'Occult blood and microbial contamination of laryngoscopes'**
M.S. Ballin et al – Stepping Hill Hospital, Stockport
5. **'Cross infection and the laryngoscope'**
David Pidduck – British Journal of Perioperative Nursing Vol 12, No5, May 2002
6. **'Infection Control in Anaesthesia'** – AAGBI Nov 2002

IS5.10 Issue 3 AAC 10.05



INTERSURGICAL®
SPECIALISTE DE L'ABORD RESPIRATOIRE

Intersurgical Sarl, 38 rue Roger Salengro - 94124 Fontenay-sous-Bois Cedex

T: 01 48 76 72 30 F: 01 48 76 91 38 info@intersurgical.fr | www.intersurgical.com

UK
info@intersurgical.com

Deutschland
info@intersurgical.de

España
info@intersurgical-es.com

Россия
info@intersurgical.ru

South Africa
info@intersurgical.co.za

Portugal
info@intersurgical.pt

Nederland
info@intersurgical.nl

Lietuva
info@intersurgical.lt

Česká Republika
info@intersurgical.cz

Philippines
info@intersurgical.ph